

NUP98 Recombinant Rabbit mAb

Catalog Number: bsm-61197R

Target Protein: NUP98

Form: Liquid

Host: Rabbit

Clonality: Recombinant

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:50-200), IHC-F (1:50-200), ICC/IF (1:50-200), IF (1:50-200)

Reactivity: Human

Predicted MW: 198

Detected MW: 188

Entrez Gene: 4928

Swiss Prot: P52948


Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: Plays a role in the nuclear pore complex (NPC) assembly and/or maintenance. Nup98 and Nup96 are involved in the bidirectional transport across the NPC. May anchor NUP153 and TPR to the NPC.

VALIDATION IMAGES



Western blot analysis of HeLa cell lysate. Using NUP98 (bsm-61197R) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.