

bs-0283R**[Primary Antibody]****SERPINB14 | Ovalbumin Rabbit pAb****BioSS**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:5000-20000)
Clonality: Polyclonal		Reactivity: Chicken
GeneID: 396058	SWISS: P01012	
Target: SERPINB14 Ovalbumin		
Purification: affinity purified by Protein A		Predicted MW.: 43 kDa
Concentration: 1mg/ml		Observed MW.: 43 kDa
Storage: 0.01M TBS (pH7.4) with 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Subcellular Location: Secreted

— VALIDATION IMAGES —

25 ug total protein per lane of various lysates
(see on figure) probed with Ovalbumin
polyclonal antibody, unconjugated (bs-0283R) at
1:10000 dilution and 4°C overnight incubation.
Followed by conjugated secondary antibody
incubation at r.t. for 60 min.

— SELECTED CITATIONS —

- **[IF=26.8]** Huiqun Zhou. et al. 3D Bioprinted Human Synovium-Cartilage Models Mimic Rheumatoid Arthritis Microenvironment and Recapitulate In Vivo Therapeutic Responses. ADV MATER. 2025 Dec;;e13952 IF ;. 41388589
- **[IF=4.493]** Gong Zheng. et al. Calreticulin as an Adjuvant In Vivo to Promote Dendritic Cell Maturation and Enhance Antigen-Specific T Lymphocyte Responses against Melanoma. J IMMUNOL RES. 2022;2022:8802004 WB ;Escherichia coli. 35983078