


bsm-61727R**[Primary Antibody]****Bioss**
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

FHL1 Recombinant Rabbit mAb**DATASHEET****Host:** Rabbit**Isotype:** IgG**Clonality:** Recombinant**CloneNo.:** 6F1**GeneID:** 2273**SWISS:** Q13642**Target:** FHL1**Immunogen:** A synthesized peptide derived from human FHL 1: 8-52.**Purification:** affinity purified by Protein A**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.**Background:** May have an involvement in muscle development or hypertrophy.**Applications:** **WB** (1:500-2000)
IP (1:20-50)**Reactivity:** Human (predicted: Mouse, Rat)**Predicted MW.:** 36**Subcellular Location:** Cytoplasm ,Nucleus**VALIDATION IMAGES**

25 ug total protein per lane of various lysates (see on figure) probed with FHL1 monoclonal antibody, unconjugated (bsm-61727R) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.