

PSMD4 Recombinant Rabbit mAb

Catalog Number: bsm-61859R

Target Protein: PSMD4

Form: Liquid

Host: Rabbit

Clonality: Recombinant

Isotype: IgG

Applications: WB (1:500-1:2000), IHC-P (1:100-1:500), IHC-F (1:100-1:500), Flow-Cyt (1:20-1:100), IF (1:100-1:500), IP (1:20-1:50)

Reactivity: Human, Mouse, Rat

Predicted MW: 41

Detected MW: 50

Entrez Gene: 5710

Swiss Prot: P55036


Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: Binds and presumably selects ubiquitin-conjugates for destruction. Displays selectivity for longer polyubiquitin chains. Modulates intestinal fluid secretion.

VALIDATION IMAGES



Western blot analysis of K562 cell lysate. Using PSMD4 (bsm-61859R) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.