

bsm-61989R**[Primary Antibody]****Beta-casein Recombinant Rabbit mAb****BioSS**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) IP (1:20-50) Reactivity: Human Predicted MW.: 25 Subcellular Location: Secreted
Clonality: Recombinant	CloneNo.: 11B2	
GeneID: 1447	SWISS: P05814	
Target: Beta-casein		
Immunogen: A synthesized peptide derived from human beta Casein: 100-150/226.		
Purification: affinity purified by Protein A		
Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.		
Background: Important role in determination of the surface properties of the casein micelles.		

— VALIDATION IMAGES —

25 ug total protein per lane of various lysates (see on figure) probed with Beta-casein monoclonal antibody, unconjugated (bsm-61989R) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.