

[Primary Antibody]

RAP1A/1B Recombinant Rabbit mAb



BioSS
ANTIBODIES

www.bioss.com.cn

sales@bioass.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET

Host: Rabbit

Isotype: IgG

Clonality: Recombinant

GeneID: 5908

SWISS: P61224,P62834

Target: RAP1A/1B

Immunogen: A synthesized peptide derived from human RAP1A: 100-150/184.

Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: RAP1A and RAP1B belong to a family of RAS-related proteins. These proteins share approximately 50% amino acid identity with the classical RAS proteins and have numerous structural features in common.

Applications: WB (1:500-2000)


Flow-Cyt (1:50-100)

IP (1:20-50)

Reactivity: Human, Mouse, Rat

Predicted MW.: 21

— VALIDATION IMAGES



25 ug total protein per lane of various lysates (see on figure) probed with RAP1A/1B monoclonal antibody, unconjugated (bsm-62137R) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.