

CD155 Recombinant Rabbit mAb

Catalog Number: bsm-62186R

Target Protein: CD155

Form: Liquid

Host: Rabbit

Clonality: Recombinant

Isotype: IgG

Applications: WB (1:500-2000), IP (1:20-50)

Reactivity: Human

Predicted MW: 45

Detected MW: 70

Entrez Gene: 5817

Swiss Prot: P15151


Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: Mediates NK cell adhesion and triggers NK cell effector functions. Binds two different NK cell receptors: CD96 and CD226. These interactions accumulates at the cell-cell contact site, leading to the formation of a mature immunological synapse between NK cell and target cell.

VALIDATION IMAGES



25 ug total protein per lane of various lysates (see on figure) probed with CD155 monoclonal antibody, unconjugated (bsm-62186R) at 1:2000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.