
PGAM1 Recombinant Rabbit mAb

Catalog Number: bsm-60524R

Target Protein: PGAM1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Recombinant

Clone No.: D3C4

Isotype: IgG/Kappa

Applications: WB (1:500-1:1000), ICC/IF (1:50)

Reactivity: (predicted:Human,Rat)

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: PGAM1 belongs to the phosphoglycerate mutase family, BPG-dependent PGAM subfamily. It is involved in the interconversion of 3- and 2-phosphoglycerate with 2,3-bisphosphoglycerate as the primer of the reaction.