

NUP62 Recombinant Rabbit mAb

Catalog Number: bsm-62402R

Target Protein: NUP62

Form: Liquid

Host: Rabbit

Clonality: Recombinant

Clone No.: 9F19

Isotype: IgG

Applications: WB (1:500-1:2000)

Reactivity: Human

Predicted MW: 53

Detected MW: 70

Entrez Gene: 23636

Swiss Prot: P37198

Purification: affinity purified by Protein A


Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: The nuclear pore complex is a structure that extends across the nuclear envelope and regulates the flow of macromolecules between the cytoplasm and the nucleus.

Nucleoporins are the main components of the nuclear pore complex in eukaryotic cells.

VALIDATION IMAGES



25 ug total protein per lane of various lysates (see on figure) probed with NUP62 monoclonal antibody, unconjugated (bsm-62402R) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.