

ACAT2 Recombinant Rabbit mAb

Catalog Number: bsm-62494R

Target Protein: ACAT2

Form: Liquid

Host: Rabbit

Clonality: Recombinant

Isotype: IgG

Applications: WB (1:1000-1:2000), IHC-P (1:100-1:200), IHC-F (1:100-1:200), IF (1:100-1:200)

Reactivity: Human, Mouse, Rat

Predicted MW: 41

Detected MW: 38

Entrez Gene: 39

Swiss Prot: Q9BWD1

Source: A synthesized peptide derived from human ACAT2: 200-397.


Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: Involved in the biosynthetic pathway of cholesterol.

VALIDATION IMAGES



Western blot analysis of K562 cell lysate. Using ACAT2 (bsm-62494R) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.