


bsm-33007R**[Primary Antibody]****BioSS**
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

GST tag Recombinant Rabbit mAb**— DATASHEET —****Host:** Rabbit**Clonality:** Recombinant**CloneNo.:** 1C14**Target:** GST tag**Purification:** affinity purified by Protein A**Storage:** 0.01M TBS(pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.**Background:** GST is a 26kDa protein encoded by the parasitic helminth *Schistosoma japonicum* and widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins.**Applications:** **WB** (1:5000-100000)**IHC-P** (1:200-1000)**IHC-F** (1:200-1000)**IF** (1:200-1000)**Flow-Cyt** (1ug/Test)**ICC/IF** (1:100-500)**ELISA** (1:1000-5000)**Reactivity:** Species independent**Observed
MW.:** 41 kDa**Subcellular
Location:** Cytoplasm**— VALIDATION IMAGES —**

100 ng rms100A7-GST protein (bs-42023P) per lane probed with GST tag monoclonal antibody respectively, unconjugated (bsm-33007R) at 1:50000 dilution and 4°C overnight incubation. Followed by corresponding conjugated secondary antibody incubation at r.t. for 60 min.