

**bsm-33034R****[ Primary Antibody ]**
**BioSS**  
**ANTIBODIES**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

## Beta tubulin Recombinant Rabbit mAb, Loading Control

### DATASHEET

**Host:** Rabbit**Clonality:** Recombinant**GeneID:** 203068**Target:** Beta tubulin**Purification:** affinity purified by Protein A**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Background:** This gene encodes a beta tubulin protein. This protein forms a dimer with alpha tubulin and acts as a structural component of microtubules. Mutations in this gene cause cortical dysplasia, complex, with other brain malformations 6. Alternative splicing results in multiple splice variants. There are multiple pseudogenes for this gene on chromosomes 1, 6, 7, 8, 9, and 13. [provided by RefSeq, Jun 2014]

**Isotype:** IgG**CloneNo.:** 7G3**SWISS:** P07437**Applications:** WB (1:50000-200000)**IHC-P** (1:200-1000)**IHC-F** (1:200-1000)**IF** (1:200-1000)**Flow-Cyt** (1μg/Test)**ICC/IF** (1:50-200)


**Reactivity:** Human, Mouse, Rat  
(predicted: Rabbit, Pig, Sheep, Cow, Chicken, Dog, GuineaPig, Monkey)

**Predicted MW.:** 50 kDa


**Observed MW.:** 50 kDa

**Subcellular Location:** Cytoplasm


### VALIDATION IMAGES




25 ug total protein per lane of various lysates (see on figure) probed with Beta tubulin monoclonal antibody, unconjugated (bsm-33034R) at 1:50000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.




Paraformaldehyde-fixed, paraffin embedded Human Breast Cancer; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation with Beta tubulin Monoclonal Antibody, Unconjugated (bsm-33034R) at 1:200 overnight at 4°C, followed by conjugation to the SP Kit (Rabbit, SP-0023) and DAB (C-0010) staining.




Paraformaldehyde-fixed, paraffin embedded Human Left Parietal Lobe; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation with Beta tubulin Monoclonal Antibody, Unconjugated (bsm-33034R) at 1:200 overnight at 4°C, followed by conjugation to the SP Kit (Rabbit, SP-0023) and DAB (C-0010) staining.




4% Paraformaldehyde-fixed HeLa (H) cell; Triton X-100 at r.t. for 20 min; Antibody incubation with (Beta tubulin) monoclonal Antibody, unconjugated (bsm-33034R) 1:100, 90 min at 37°C; followed by BF594 conjugated Goat Anti-Rabbit IgG antibody (red, bs-60295G-BF594) at 37°C for 90 min, then stained with Phalloidin/FITC for 30 min at room temperature; DAPI (blue, C02-04002) was used to stain the cell nuclei.



4% Paraformaldehyde-fixed HeLa (H) cell; Triton X-100 at r.t. for 20 min; Antibody incubation with (Beta tubulin) monoclonal Antibody, unconjugated (bsm-33034R) 1:100, 90 min at 37°C; followed by conjugated Goat Anti-Rabbit IgG antibody (green, bs-60295G-BF488) at 37°C for 90 min, DAPI (blue, C02-04002) was used to stain the cell nuclei. PBS instead of the primary antibody was used as the blank control.



Paraformaldehyde-fixed, paraffin embedded Human Colon; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation with Beta tubulin Polyclonal Antibody, Unconjugated (bsm-33034R) at 1:200 overnight at 4°C. Followed by conjugated Goat Anti-Rabbit IgG antibody (green, bs-0295G-BF488), DAPI (blue, C02-04002) was used to stain the cell nuclei.



The HeLa (H) cells were fixed with 4% PFA (10 min at r.t.) and then permeabilized with 90% ice-cold methanol for 20 min at -20°C, the cells then were incubated in 5% BSA to block non-specific protein-protein interactions (30 min at r.t.). Primary Antibody (green): Rabbit Anti-Beta tubulin antibody (bs-33034R): 1 µg/10<sup>6</sup> cells; Secondary Antibody (white blue): Goat anti-Rabbit IgG-BF488(bs-60295G-BF488): 1 µg/test. Isotype Control (orange): Rabbit IgG (bs-0295P). Blank control (black): PBS. Acquisition of 20,000 events was performed.

## — SELECTED CITATIONS —

- **[IF=17.1]** Lei Liu. et al. Myricetin Oligomer Triggers Multi-Receptor Mediated Penetration and Autophagic Restoration of Blood-Brain Barrier for Ischemic Stroke Treatment. ACS NANO. 2024;XXXX(XXX):XXX-XXX WB ;Rat,Mouse. 38533773
- **[IF=11.9]** Xinying Wang. et al. An Exosome-Inspired Docetaxel Prodrug Nanoplatfrom for Potent STING Activation and Synergistic Chemoimmunotherapy. ASIAN J PHARM SCI. 2026 Feb;;101139 WB ;MOUSE. 10.1016/j.ajps.2026.101139