
Mouse Anti-Rabbit IgM, Cy3 conjugated

Catalog Number: bs-0369M-Cy3

Target Protein: Mouse Anti-Rabbit IgM

Concentration: 2.0 mg/ml

Form: Liquid

Host: Mouse

Clonality: Polyclonal

Isotype: IgG

Applications: IF (1:100-1000)

Excitation spectrum: 514nm,552nm

Emission spectrum: 570nm

Reactivity: Rabbit

Purification: affinity purified by Protein G

Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol.

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Immunoglobulin M (IgM) normally constitutes about 10% of serum immunoglobulins. IgM antibody is prominent in early immune responses to most antigens and is largely confined to plasma due to its large size. Monomeric IgM is expressed as a membrane bound antibody on the surface of B cells and as a pentamer when secreted by plasma cells. Due to its high valency IgM is more efficient than other isotypes in binding antigens with repeating epitopes (virus particles and red blood cells) and is more efficient than IgG in activating the complement pathway. The gene for the mu constant region contains four domains separated by short intervening sequences.

PRODUCT SPECIFIC PUBLICATIONS

[IF=5.811] Lingli Ding, et al. Ginsenoside Compound K Enhances Fracture Healing via Promoting Osteogenesis and Angiogenesis.. FRONT PHARMACOL. 2022 Apr;13:855393-855393 IF ; Rat . 35462912

[IF=4.7] Ding Lingli, et al. Ginsenoside compound-K attenuates OVX-induced osteoporosis via the suppression of RANKL-induced osteoclastogenesis and oxidative stress. NAT PRODUCT BIOPROSP. 2023 Dec;13(1):1-12 IF ; Mouse . 37940733