

bs-18471R

[Primary Antibody]

ZDHHC1/ZNF377 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000)
Clonality: Polyclonal		IHC-P (1:100-500)
GeneID: 29800	SWISS: Q8WTX9	IHC-F (1:100-500)
Target: ZDHHC1/ZNF377		IF (1:100-500)
Immunogen: KLH conjugated synthetic peptide derived from human ZDHHC1/ZNF377: 1-100/485.		ICC/IF (1:100-500)
Purification: affinity purified by Protein A		Reactivity: (predicted: Human, Mouse, Rat)
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Predicted MW.: 55 kDa
		Subcellular Location: Cell membrane