

bs-0032R**[Primary Antibody]****Bioss**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn


400-901-9800

Bcl-2 Rabbit pAb**— DATASHEET —**


Host: Rabbit	Isotype: IgG	Applications: WB (1:1000-5000) IHC-P (1:200-500) IHC-F (1:200-500) IF (1:200-500) ELISA (1:5000-10000) Reactivity: Human, Mouse, Rat (predicted: Chicken) Predicted MW.: 26 kDa Observed MW.: 26 kDa Subcellular Location: Cell membrane ,Cytoplasm ,Nucleus
Clonality: Polyclonal		
GeneID: 596	SWISS: P10415	
Target: Bcl-2		
Immunogen: KLH conjugated synthetic peptide derived from human Bcl-2: 101-160/236.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: The Bcl-2 gene was isolated at the chromosomal breakpoint of t(14;18)-bearing follicular B cell lymphomas(1,2).Bcl-2 blocks cell death following a variety of stimuli and confers a death-sparing effect to certain hematopoietic cell lines following growth factor withdrawal (3,5).Bcl-2 appears to function in several subcellular locations yet lacks any known motifs that would confer insight into its mechanism of action (6,7).A more recently identified protein,designated Bax p21(i.e., Bcl-associated X protein),has extensive amino acid homology with Bcl-2 and both homodimerizes and forms heterodimers with Bcl-2(8). Overexpression of Bax accelerates apoptotic death induced by cytokine deprivation in an IL-3 dependent cell line and Bax also counters the death repressor activity of Bcl-2(8).		

— VALIDATION IMAGES —

25 ug total protein per lane of various lysates (see on figure) probed with Bcl-2 polyclonal antibody, unconjugated (bs-0032R) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.



Paraformaldehyde-fixed, paraffin embedded Human Breast Cancer; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with Bcl-2 Polyclonal Antibody, Unconjugated (bs-0032R) at 1:200 overnight at 4°C, followed by conjugation to the bs-0295G-HRP and DAB (C-0010) staining.



Paraformaldehyde-fixed, paraffin embedded Human Lymph; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with Bcl-2 Polyclonal Antibody, Unconjugated (bs-0032R) at 1:200 overnight at 4°C, followed by conjugation to the bs-0295G-HRP and DAB (C-0010) staining.

— SELECTED CITATIONS —

- **[IF=18.027]** Shikai Liu. et al. On-Demand Generation of Peroxynitrite from an Integrated Two-Dimensional System for Enhanced Tumor Therapy. ACS NANO. 2022;XXXX(XXX):XXX-XXX WB ;Human. 35666853
- **[IF=16.6]** Liu Chunxiao. et al. Targeting P2Y14R protects against necroptosis of intestinal epithelial cells through PKA/CREB/RIPK1 axis in ulcerative colitis. NAT COMMUN. 2024 Mar;15(1):1-16 WB ;Mouse. 38453952
- **[IF=14.7]** Chunxiao Liu. et al.Targeting P2Y14R protects against necroptosis of intestinal epithelial cells through PKA/CREB/RIPK1 axis in ulcerative colitis.nature communications.2024 Mar 7;15(1):2083. Western blot ;Mouse.

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

38453952

- **[IF=14.593]** Guo-Bin Ding. et al. Molecularly engineered tumor acidity-responsive plant toxin gelonin for safe and efficient cancer therapy. BIOACT MATER. 2022 Dec;18:42 IHC ;Mouse. 10.1016/j.bioactmat.2022.02.001
- **[IF=14.3]** Huitong Luo. et al. Nanocomposite hydrogel orchestrating multiple modulation of degenerative microenvironment for potential application in intervertebral disc regeneration. J MATER SCI TECHNOL. 2025 Feb;; WB,IP ;Human. 10.1016/j.jmst.2024.12.042