

bs-41382R

[Primary Antibody]

ASFV p30 Rabbit pAb

BioSS
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn


techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ELISA (1:5000-10000)
Clonality: Polyclonal		Reactivity: ASFV
Target: ASFV p30		
Immunogen: Recombinant ASFV p30 protein: 1-204/204.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		Predicted MW.: 23.5 kDa
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Observed MW.: 32 kDa

— VALIDATION IMAGES —



Sample: Lane 1: Recombinant ASFV p30 protein,
N-His Primary: Anti- ASFV p30 (bs-41382R) at
1/1000 dilution Secondary: IRDye800CW Goat
Anti-Rabbit IgG at 1/20000 dilution Predicted
band size: 23.5 kDa Observed band size: 32 kDa

— SELECTED CITATIONS —

- **[IF=6.1]** Yuhao Huang. et al. Screening of aptamer for p30 protein of African swine fever virus and establishment of electrochemical immunosensor for detection of p30 antibody. TALANTA. 2025 Mar;;127964 ;. 40154044