

Gastrin 17 (G-17) Mouse mAb

Cat. No:	V4001
Name:	Gastrin 17 (G-17) Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	9F8
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	Human
Format:	Liquid
Concentration:	≥1mg/ml
Buffer:	0.01M PBS (pH7.4) with 0.02% Proclin300.
Applications:	CLIA
Antigen:	Recombinant human Gastrin 17
Pairing Information:	

CL		Detection	
		V4001	V4002
Capture	V4001	-	+++
	V4002	-	-