

bsm-52727R**[Primary Antibody]****BioSS**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

IL-1 beta Recombinant Rabbit mAb**— DATASHEET —**

Host: Rabbit Clonality: Recombinant GeneID: 3553 Target: IL-1 beta Immunogen: A synthesized peptide derived from human IL 1 beta: 200-269. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 0.05% BSA, 40% Glycerol, 0.02% Proclin300. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: The protein encoded by this gene is a member of the interleukin 1 cytokine family. This cytokine is produced by activated macrophages as a proprotein, which is proteolytically processed to its active form by caspase 1 (CASP1/ICE). This cytokine is an important mediator of the inflammatory response, and is involved in a variety of cellular activities, including cell proliferation, differentiation, and apoptosis. The induction of cyclooxygenase-2 (PTGS2/COX2) by this cytokine in the central nervous system (CNS) is found to contribute to inflammatory pain hypersensitivity. This gene and eight other interleukin 1 family genes form a cytokine gene cluster on chromosome 2. [provided by RefSeq, Jul 2008].	Isotype: IgG CloneNo.: 1T3 SWISS: P01584	Applications: WB (1:1000-5000) Reactivity: Human Predicted MW.: 31 kDa Observed MW.: 35 kDa Subcellular Location: Secreted
--	---	---

— VALIDATION IMAGES —

THP-1 (H) cells were treated with or without TPA (80nM) for 24h, then treated with (100ng/ml) lipopolysaccharide (LPS) for 6h, then with (300ng/ml) Brefeldin A added after 3h, 25 µg total protein per lane of cell lysates (see on figure) probed with IL-1 beta monoclonal antibody, unconjugated (bsm-52727R) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.

— SELECTED CITATIONS —

- **[IF=8.2]** Lin Peng. et al. Bisphenol A exposure exacerbates tracheal inflammatory injury in selenium-deficient chickens by regulating the miR-155/TRAFF3/ROS pathway. INT J BIOL MACROMOL. 2023 Dec;253:127501 WB ;Chicken. 37866585
- **[IF=5.7]** Mengyue Yin. et al. The impact of SARS-CoV-2 infection at different stages of pregnancy on placental inflammatory responses. Frontiers in Immunology. Western blot ;Human. 33613576
- **[IF=5.1]** Min Wei. et al. The role of ROS-pyoptosis in PM2.5 induced air-blood barrier destruction. CHEM-BIOL INTERACT. 2023 Dec;386:110782 WB ;Human. 37884181

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

- **[IF=2.6]** Yu Fang. et al. Corydalis yanhusuo Polysaccharides regulates HPA-axis mediated microglia activation and inhibits astrocyte A1 transformation to improve depression-like behavior. BRAIN RES. 2025 Oct;1864:149780 WB ;Mouse. 40513699
- **[IF=preprint]** Mengyue Yin. et al.Impact of SARS-CoV-2 Infection During Various Pregnancy Trimesters on Maternal and Fetal Outcomes.Research Square. Western blot ;Human. 10.21203/rs.3.rs-4488503/v1