


bsm-52384R**[Primary Antibody]****CD63 Recombinant Rabbit mAb****BioSS**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn


400-901-9800

— DATASHEET —**Host:** Rabbit**Isotype:** IgG**Clonality:** Recombinant**CloneNo.:** 4C2**GeneID:** 967**SWISS:** P08962**Target:** CD63**Immunogen:** A synthesized peptide derived from human CD63: 100-200.**Purification:** affinity purified by Protein A**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.**Background:** The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. The encoded protein is a cell surface glycoprotein that is known to complex with integrins. It may function as a blood platelet activation marker. Deficiency of this protein is associated with Hermansky-Pudlak syndrome. Also this gene has been associated with tumor progression. Alternative splicing results in multiple transcript variants encoding different protein isoforms. [provided by RefSeq, Apr 2012]**Applications:** WB (1:1000-5000)**IHC-P** (1:200-800)**IHC-F** (1:200-800)**IF** (1:200-800)**ICC/IF** (1:50-200)**Reactivity:** Human**Predicted MW.:** 26 kDa**Observed MW.:** 30-60 kDa**Subcellular Location:** Cell membrane ,Cytoplasm**— VALIDATION IMAGES —**


25 ug total protein per lane of various lysates (see on figure) probed with CD63 monoclonal antibody, unconjugated (bsm-52384R) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.



Paraformaldehyde-fixed, paraffin embedded Human Tonsil; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with CD63 Monoclonal Antibody, Unconjugated (bsm-52384R) at 1:200 overnight at 4°C, followed by conjugation to the bs-0295G-HRP and DAB (C-0010) staining.



Paraformaldehyde-fixed, paraffin embedded Human Lung Cancer; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with CD63 Monoclonal Antibody, Unconjugated (bsm-52384R) at 1:200 overnight at 4°C, followed by conjugation to the bs-0295G-HRP and DAB (C-0010) staining.



4% Paraformaldehyde-fixed A375 (H) cell; Triton X-100 at r.t. for 20 min; Antibody incubation with (CD63) monoclonal Antibody, unconjugated (bsm-52384R) 1:100, 90 min at 37°C; followed by conjugated Goat Anti-Rabbit IgG antibody (green, bs-40295G-FITC) at 37°C for 90 min, DAPI (blue, C02-04002) was used to stain the cell nuclei. PBS instead of the primary antibody was used as the blank control.

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

— SELECTED CITATIONS —

- **[IF=4.6]** Huilong Li. et al. Exosomes derived from syncytia induced by SARS-2-S promote the proliferation and metastasis of hepatocellular carcinoma cells. FRONTIERS IN CELLULAR AND INFECTION MICROBIOLOGY. 2025 Jan 8;14:1415356. Western blot ;Human. 39844837
- **[IF=4]** Liang Min. et al. Astragaloside IV Suppresses the Effects of Hepatocellular Carcinoma Cells on Proliferation, Angiogenesis, and Invasion in Human Umbilical Vein Endothelial Cells by Controlling Exosomes by Inhibiting Rab27a. J FOOD BIOCHEM. 2023;2023:8812742 WB ;Human. 10.1155/2023/8812742
- **[IF=3.9]** Irene Chavarría-Cubel. et al. The synthetic TRPML1 agonist ML-SA1 mitigates intracellular lipid accumulation induced by antipsychotics in vitro by stimulating release of extracellular microvesicles. BBA-MOL CELL BIOL L. 2025 May;1870:159611 IF ;Human. 40222411
- **[IF=3.251]** Qifang Niu. et al. Exosomes Derived from Bone Marrow Mesenchymal Stem Cells Alleviate Ischemia-Reperfusion Injury and Promote Survival of Skin Flaps in Rats. LIFE-BASEL. 2022 Oct;12(10):1567 WB ;Rat. 36295004
- **[IF=3.21]** Ming-Zhi Huang. et al. Exosomes from artesunate-treated bone marrow-derived mesenchymal stem cells transferring SNHG7 to promote osteogenesis via TAF15-RUNX2 pathway. REGEN MED. 2022 Oct 02 WB ;Mouse, Human. 36184881