

bsm-52283R**[Primary Antibody]****BioSS**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn


400-901-9800

SCARB1/Scavenger Receptor BI Recombinant Rabbit mAb**DATASHEET**


Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:200-500) IF (1:100-500) ICC/IF (1:50-200) Reactivity: Human, Mouse, Rat Predicted MW.: 61 kDa Observed MW.: 75 kDa Subcellular Location: Cell membrane ,Cytoplasm
Clonality: Recombinant	CloneNo.: 1C2	
GeneID: 949	SWISS: Q14108	
Target: SCARB1/Scavenger Receptor BI		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: High density lipoproteins (HDLs) play a critical role in cholesterol metabolism and their plasma concentrations are inversely correlated with risk for atherosclerosis. The SR-BI (Scavenger Receptor BI) protein binds HDLs and mediates selective uptake of HDL cholesteryl ester. SR-BI binds HDL with high affinity, is expressed primarily in liver and nonplacental steroidogenic tissues, and mediates selective cholesterol uptake by a distinct mechanism. In mice, it seems that SR-BI plays a key role in determining the levels of plasma lipoprotein cholesterol and the accumulation of cholesterol stores in the adrenal gland. Scavenging Receptor SR-BI plays a critical role in HCV attachment and/or cell entry by interacting with HCV E1/E2 glycoproteins heterodimer.		

VALIDATION IMAGES


Sample: Lane 1: Mouse Liver tissue lysates Lane 2: Mouse Adrenal gland tissue lysates Lane 3: Rat Liver tissue lysates Lane 4: Rat Adrenal gland tissue lysates Lane 5: Human HepG2 cell lysates Lane 6: Human HeLa cell lysates Primary: Anti-SCARB1/Scavenger Receptor BI (bsm-52283R) at 1/500 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 61 kDa Observed band size: 75 kDa



Paraformaldehyde-fixed, paraffin embedded (mouse liver); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Incubation with (SCARB1/Scavenger Receptor BI) Monoclonal Antibody, Unconjugated (bsm-52283R) at 1:2000 overnight at 4°C, followed by operating according to SP Kit(Rabbit) (sp-0023) instructions and DAB staining.



Paraformaldehyde-fixed, paraffin embedded (mouse adrenal gland); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Incubation with (SCARB1/Scavenger Receptor BI) Monoclonal Antibody, Unconjugated (bsm-52283R) at 1:2000 overnight at 4°C, followed by operating according to SP Kit(Rabbit) (sp-0023) instructions and DAB staining.



Cell line: HepG2 Fixation: 4% Paraformaldehyde Permeabilization: 0.1% TritonX-100 Primary Ab dilution: 1:50 Primary Ab incubation condition:

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

4°C overnight Secondary Ab: Goat Anti-Rabbit
IgG Nuclear counter stain: DAPI (Blue) Comment:
Color green is the positive signal for bsm-52283R

— SELECTED CITATIONS —

- **[IF=5.8]** Min-Chien Tsai. et al. Cav3.1 T-type calcium channel blocker NNC 55-0396 reduces atherosclerosis by increasing cholesterol efflux. BIOCHEM PHARMACOL. 2024 Apr;222:116096 WB ;Human. 38423188
- **[IF=3.9]** Zhou Qianhui. et al. SEC14L2 regulates the transport of cholesterol in non-small cell lung cancer through SCARB1. LIPIDS HEALTH DIS. 2024 Dec;23(1):1-13 IF ;Human. 39696431