
Rabbit Anti-Dog IgM

Catalog Number: bs-0373R

Target Protein: Rabbit Anti-Dog IgM

Form: Lyophilized

Host: Rabbit

Clonality: Polyclonal

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:50-200), ELISA (1:500-1000)

Reactivity: Dog (predicted:)

Isotype: IgG

Source: Native Dog IgM: full length.

Purification: affinity purified by Protein A

Storage: 0.01M PBS (pH7.4)

Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.