

Recombinant PRRSV N protein, N-His

Catalog Number: bs-41387P

AA Seq: 1-128/128

Predicted MW: 13.7

Detected MW: 18 kDa

Tags: N-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH9.0) with 8M urea.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: PRRSV is a small, enveloped RNA virus. It contains a single-stranded, positive-sense, RNA genome with a size of approximately 15 kilobases. The genome contains nine open reading frames. PRRSV is a member of the genus Arterivirus, family Arteriviridae, order Nidovirales. Subclinical infections are common, with clinical signs occurring sporadically in a herd. Clinical signs include reproductive failure in sows such as abortions and giving birth to stillborn or mummified fetuses, and cyanosis of the ear and vulva. In neonatal pigs, the disease causes respiratory distress, with increased susceptibility to respiratory infections such as Glasser's disease.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.