

bs-0295M-BF594

[Secondary Antibodies]

Mouse Anti-Rabbit IgG H&L, BF594 conjugated

Bioss
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Mouse	Isotype: IgG	Applications: IF (1:100-1000)
Clonality: Polyclonal		Flow-Cyt (1:100-1000)
Target: Mouse Anti-Rabbit IgG H&L		ICC/IF (1:100-1000)
Purification: affinity purified by Protein G		Excitation Spectrum: 594nm
Concentration: 2.0 mg/ml		Emission spectrum: 621nm
Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Reactivity: Rabbit

— SELECTED CITATIONS —

- **[IF=4.9]** Qian Zhang. et al. Proteomics and Expression of HIF2 α /BNIP3L Signaling in Yak Brains at Different Altitudes. INTERNATIONAL JOURNAL OF MOLECULAR SCIENCES. 2025 Feb 16;26(4):1675. IF ;Yak. 10.3390/ijms26041675
- **[IF=4.9]** Qian Zhang. et al. Proteomics and Expression of HIF2 α /BNIP3L Signaling in Yak Brains at Different Altitudes. int. j. mol. sci.. 2025 Feb 16;26(4):1675. IF ;yak. 40004139