

bsm-60750R**[Primary Antibody]****BioSS**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn


400-901-9800

TBK1 Recombinant Rabbit mAb**DATASHEET**


Host: Rabbit	Isotype: IgG	Applications: WB (1:1000-5000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) Reactivity: Human, Mouse, Rat Predicted MW.: 84 kDa Observed MW.: 84 kDa Subcellular Location: Cytoplasm
Clonality: Recombinant	CloneNo.: B5H9	
GeneID: 29110	SWISS: Q9UHD2	
Target: TBK1		
Immunogen: A synthesized peptide derived from human TBK1: 1-39.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: The NF-kappa-B (NFKB) complex of proteins is inhibited by I-kappa-B (IKB) proteins, which inactivate NFKB by trapping it in the cytoplasm. Phosphorylation of serine residues on the IKB proteins by IKB kinases marks them for destruction via the ubiquitination pathway, thereby allowing activation and nuclear translocation of the NFKB complex. The protein encoded by this gene is similar to IKB kinases and can mediate NFKB activation in response to certain growth factors. [provided by RefSeq, Oct 2010]		

VALIDATION IMAGES

25 ug total protein per lane of various lysates (see on figure) probed with TBK1 monoclonal antibody, unconjugated (bsm-60750R) at 1:2000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.



Paraformaldehyde-fixed, paraffin embedded Human Pancreatic Cancer; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation with TBK1 Monoclonal Antibody, Unconjugated(bsm-60750R) at 1:200 overnight at 4°C, followed by conjugation to the SP Kit (Rabbit, SP-0023) and DAB (C-0010) staining.



Paraformaldehyde-fixed, paraffin embedded Human Placenta; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation with TBK1 Monoclonal Antibody, Unconjugated(bsm-60750R) at 1:200 overnight at 4°C, followed by conjugation to the SP Kit (Rabbit, SP-0023) and DAB (C-0010) staining.

SELECTED CITATIONS

- **[IF=4.8]** Takumi Ikeda. et al. Molecular Cloning, Expression Analyses, and Physiological Roles of Cathelicidins in the Bursa of Fabricius of the Japanese Quail, Coturnix japonica. ANTIBIOTICS-BASEL. 2023 Aug;12(8):1341 IHC ;Quail. 37627761