

Recombinant MPXV M1R protein, N-His

Catalog Number: bs-42062P

Concentration: >0.5 mg/ml

AA Seq: 1-182/250

Predicted MW: 25.1

Tags: N-His

QC Testing

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE


Purification: AC

Form: Lyophilized or Liquid


Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Monkeypox is a rare zoonosis caused by monkeypox virus, which has become the most serious orthopoxvirus and consists of complex double stranded DNA. The cases are mostly in central and western Africa. The pathogenesis of monkeypox is that the virus invades the body from respiratory mucosa, multiplies in lymphocytes, and incurs into blood producing transient venereal toxemia. After the virus multiplies in cells, the cells can invade the blood and propagate to the skin of the whole body, causing lesions.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.