

bsm-54306R**[Primary Antibody]**

Cytochrome P450 17A1 Recombinant Rabbit mAb

BioSS
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn


techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) Flow-Cyt (1:50-100) ICC/IF (1:50-200) Reactivity: Human, Mouse, Rat Predicted MW.: 57 kDa Subcellular Location: Cell membrane
Clonality: Recombinant	CloneNo.: 2F7	
GeneID: 1586	SWISS: P05093	
Target: Cytochrome P450 17A1		
Immunogen: A synthesized peptide derived from human Cytochrome P450 17A1: 101-220.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: Cytochrome P450 17A1 (CYP17A1) belongs to the cytochrome P450 family; it plays a role in the conversion of pregnenolone and progesterone into their 17-alpha-hydroxylated products and subsequently to dehydroepiandrosterone (DHEA) and androstenedione. CYP17A1 also catalyzes both the 17-alpha-hydroxylation and the 17,20-lyase reaction. CYP17A1 is involved in sexual development during fetal life and at puberty. Defects in CYP17A1 are the cause of adrenal hyperplasia type 5 (AH5). AH5 is a form of congenital adrenal hyperplasia, a common recessive disease due to defective synthesis of cortisol.		

— VALIDATION IMAGES —



25 ug total protein per lane of various lysates (see on figure) probed with Cytochrome P450 17A1 monoclonal antibody, unconjugated (bsm-54306R) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.

— SELECTED CITATIONS —

- **[IF=7.129]** Hui Zhao. et al. Alleviating effects of selenium on fluoride-induced testosterone synthesis disorder and reproduction toxicity in rats. ECOTOX ENVIRON SAFE. 2022 Dec;247:114249 WB ;Rat. 36323150
- **[IF=2.7]** Qingjing Gao. et al. Seasonal changes in endoplasmic reticulum stress and steroidogenesis in the ovary of the wild ground squirrels (Citellus dauricus Brandt). GEN COMP ENDOCR. 2023 Nov;343:114368 IHC ;Squirrel. 37604348