

bs-20521R

[Primary Antibody]

PrV gE Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: PrV gE

Immunogen: KLH conjugated synthetic peptide derived from PrV gE:
441-540/577.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%
Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated
freeze/thaw cycles.

Applications: ELISA (1:5000-10000)

Reactivity: (predicted: PrV)

Predicted MW.: 62 kDa