

Aurora B Recombinant Rabbit mAb

Catalog Number: bsm-52019R

Target Protein: Aurora B

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Recombinant

Clone No.: 1G9

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:50-200), ICC/IF (1:50-200), IF (1:50-200)

Reactivity: Human (predicted:Mouse,Rat)

Predicted MW: 39 kDa

Entrez Gene: 9212

Swiss Prot: Q96GD4


Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.


Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a member of the aurora kinase subfamily of serine/threonine kinases. The genes encoding the other two members of this subfamily are located on chromosomes 19 and 20. These kinases participate in the regulation of segregation of chromosomes during mitosis and meiosis through association with microtubules. A pseudogene of this gene is located on chromosome 8. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012].

VALIDATION IMAGES



Western blot analysis of Aurora B on different lysates using anti-Aurora B antibody at 1/1,000 dilution.
Positive control: Lane 1: Hela Lane 2: MCF-7



Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-Aurora B antibody. Counter stained with hematoxylin.

PRODUCT SPECIFIC PUBLICATIONS

[IF=1.7] Zhaofeng Zhao. et al. Role of FOXM1 and AURKB in regulating keratinocyte function in psoriasis. OPEN MED-WARSAW. 2024 Jan;19(1): IF ; Human . 39381423