

bs-0297G-Gold

[Secondary Antibodies]

Goat Anti-Human IgG H&L, Gold conjugated

BioSS
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Goat

Isotype: IgG

Clonality: Polyclonal

Target: Goat Anti-Human IgG H&L

Purification: affinity purified by Protein G, nonspecific adsorbed

Concentration: 0.4 mg/ml

Storage: 20mM TBS (pH8.0) with 1% BSA and 0.02% Proclin300.
Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.

Background: Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.

Applications: **ICA** (1:20-200)
IEM (1:20-200)

Reactivity: Human

— SELECTED CITATIONS —

- **[IF=5.2]** Feifan Zheng, et al. Development of single molecule binding kinetics immunosensor based on surface plasmon resonance microscopy. MEASUREMENT. 2025 Jan;242:116316 ;. 10.1016/j.measurement.2024.116316