


1×PBS缓冲液（干粉，2L） PBS (1×, powder, 2L)

产品编号：C0103

保存条件：室温保存.

产品介绍：PBS缓冲液（phosphate buffered solution）即磷酸盐缓冲液，能够提供相对稳定的离子环境和pH缓冲能力，是生物学中经常使用的缓冲盐溶液。pH为 7.2-7.4，主要成分为磷酸二氢钠、磷酸氢二钠以及氯化钠，常用于分子克隆及免疫组化等实验。本产品为配制2L 1×PBS缓冲液的粉剂，使用前加2L蒸馏水或去离子水溶解即可。

产品图片



PBS (1×, powder, 2L) (C0103)

PRODUCT SPECIFIC PUBLICATIONS

[IF=27.7] Yu-Zhang Wang. et al. Targeted tumor therapy with L-cyst(e)ine-addicted bacteria-nanodrug biohybrids. CELL METAB. 2025 Apr Other ; . 40215982

[IF=11.9] Zhanliang Jiang. et al. Hypoxia, Cuproptosis, and Osteoarthritis: Unraveling the Molecular Crosstalk. REDOX BIOL. 2025 Jul;:103757 ; . 40669206

[IF=9.1] Ziling Ye. et al. Molecular Imprinting-Based Thickness-Encoded Suspension Arrays. ACS SENSORS. 2025;XXXX(XXX):XXX-XXX ; . 41389041

[IF=6.656] Zhenwei Zhou. et al. BuShen JianGu Fang alleviates cartilage degeneration via regulating multiple genes and signaling pathways to activate NF-κB/Sox9 axis. PHYTOMEDICINE. 2023 May;113:154742 Other ; . 36893673

[IF=6.388] Xin Liu. et al. Astragaloside IV ameliorates radiation-induced nerve cell

damage by activating the BDNF/TrkB signaling pathway. PHYTOTHER RES. 2023

May;; Other ; . 37226643