
ELISA标准品/样品稀释液 Sample dilution buffer

产品编号: C5212

保存条件: Store at 2-8°C.

产品介绍: 有利于保持样品蛋白质的稳定性; 提高信噪比; 提高免疫实验的敏感度; 降低实验背景。

PRODUCT SPECIFIC PUBLICATIONS

[IF=16] Yaqi Huang. et al. Single-Protein Determinations by Magnetofluorescent Qubit Imaging with Artificial-Intelligence Augmentation at the Point-Of-Care. ACS NANO. 2025;XXXX(XXX):XXX-XXX ; . 40388114