

bs-0296G-SA

[Secondary Antibodies]

Goat Anti-Mouse IgG H&L, SA conjugated

Bioss
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Goat Clonality: Polyclonal Target: Goat Anti-Mouse IgG H&L Purification: affinity purified by Protein G Concentration: 1.0 mg/ml Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.	Isotype: IgG	Reactivity: Mouse
--	---------------------	--------------------------

— SELECTED CITATIONS —

- **[IF=3.4]** Jianfei Wu. et al. A potential strategy for improving offspring behavior in maternal immune activation: Amantadine-mediated suppression of neuroinflammation. Animal Models and Experimental Medicine. 2025 Jul; WB ;Mouse. 40697163