

bsm-41049M

[Primary Antibody]

FMDV VP1 protein (type O) Mouse mAb

Bioss
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn


techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Mouse	Isotype: IgG1	Applications: WB (1:500-2000) ELISA (1:5000-10000) Reactivity: FMDV
Clonality: Monoclonal	CloneNo.: 2F11	
Target: FMDV VP1 protein (type O)		
Purification: affinity purified by Protein G		
Concentration: 1mg/ml		
Storage: Size : 50ul/100ul/200ul 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Size : 200ug (PBS only) 0.01M PBS Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— VALIDATION IMAGES —



Sample: Lane 1: Recombinant FMDV VP1 (type O) protein, Sumo & His
Primary: Anti-FMDV VP1 protein (type O) (bsm-41049M) at 1/1000 dilution
Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution
Predicted band size: kDa
Observed band size: 46 kDa